

# FPGT antibody



Catalog Number: 110718

## Product name

FPGT antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FPGT Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue. Positive IHC detected in human colon tissue, human kidney tissue. Observed molecular weight by Western-blot: 38 kDa

## Preparation

This antibody was obtained by immunization of FPGT recombinant protein (Accession Number: NM\_003838). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:110718(FPGT Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:110718(FPGT Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:110718(FPGT Antibody) at dilution of 1:50

# FPGT antibody



Catalog Number: 110718

---

(under 40x lens)