

FMN2 antibody



Catalog Number: 110726

Product name

FMN2 antibody

Immunogen

Specificity

Human, Mouse; other species not tested.

Antibody description

FMN2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human gliomas tissue, human ovary tumor tissue. Observed molecular weight by Western-blot: 195 kDa

Preparation

This antibody was obtained by immunization of FMN2 recombinant protein (Accession Number: BC014364). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110726(FMN2 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110726(FMN2 antibody) at dilution of 1:100 (under 40x lens)