

FN3K antibody

Catalog Number: 110734

Product name

FN3K antibody

WB: 1:1000-1:10000

IP: 1:200-1:2000

Specificity


Human, Mouse; other species not tested.

IF: 1:50-1:500

Antibody description

FN3K Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue, human heart tissue, human kidney tissue, mouse skeletal muscle tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 35 kDa

Validations



Preparation

This antibody was obtained by immunization of FN3K recombinant protein (Accession Number: NM_022158). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110734(FN3K antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IP, IF

Dilutions

Recommended Dilution:



IP Result of anti-FN3K (IP:Catalog No:110734, 3ug; Detection:Catalog No:110734 1:500) with HeLa cells lysate 1840ug.



FN3K antibody



Catalog Number: 110734

Immunofluorescent analysis of HeLa cells using
Catalog No:110734(FN3K Antibody) at dilution of

1:50 and Alexa Fluor 488-conjugated AffiniPure
Goat Anti-Rabbit IgG(H+L)