

# FXYD4 antibody

Catalog Number: 110756

## Product name

FXYD4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FXYD4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, human kidney tissue, human testis tissue, mouse kidney tissue. Positive IHC detected in mouse testis tissue. Observed molecular weight by Western-blot: 9-10 kDa

## Preparation

This antibody was obtained by immunization of FXYD4 recombinant protein (Accession Number: NM\_001184963). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110756(FXYD4 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog No:110756(FXYD4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog

# **FXVD4 antibody**



Catalog Number: 110756

---

No:110756(FXVD4 Antibody) at dilution of 1:50 (under 40x lens)