

G6PD antibody

Catalog Number: 110766

Product name

G6PD antibody

Specificity

Human; other species not tested.

Antibody description

G6PD Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 59 kDa

Preparation

This antibody was obtained by immunization of G6PD recombinant protein (Accession Number: NM_001360016). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110766(G6PD Antibody) at dilution of 1:1000