

FTSJ2 antibody

Catalog Number: 110787

Product name

FTSJ2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FTSJ2 Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in human skeletal muscle tissue, human placenta tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of FTSJ2 recombinant protein (Accession Number: NM_013393). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110787(FTSJ2 Antibody) at dilution of 1:600



IP Result of anti-FTSJ2 (IP:Catalog No:110787, 3ug; Detection:Catalog No:110787 1:500) with mouse skeletal muscle tissue lysate 4000ug.