

# GADD45A antibody



Catalog Number: 110820

## Product name

GADD45A antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GADD45A Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in PC-13 cells, A549 cells, COLO 320 cells, K-562 cells, mouse lung tissue, mouse spleen tissue, rat lung tissue. Observed molecular weight by Western-blot: 18 kDa or 36 kDa

## Preparation

This antibody was obtained by immunization of GADD45A recombinant protein (Accession Number: NM\_001924). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



PC-13 cells were subjected to SDS PAGE followed by western blot with Catalog No:110820(GADD45A Antibody) at dilution of



1:600

IP Result of anti-GADD45A (IP:Catalog No:110820, 4ug; Detection:Catalog No:110820 1:500) with K-562 cells lysate 3200ug.