

GAL3ST4 antibody



Catalog Number: 110830

Product name

GAL3ST4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GAL3ST4 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, A2780 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 55-60kd

Preparation

This antibody was obtained by immunization of GAL3ST4 recombinant protein (Accession Number: NM_024637). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110830(GAL3ST4 antibody) at dilution of 1:500