

GATM antibody



Catalog Number: 110887

Product name

GATM antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GATM Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse kidney tissue, human kidney tissue, rat kidney tissue. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of GATM recombinant protein (Accession Number: NM_001482). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110887(GATM antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110887(GATM antibody) at dilution of 1:100 (under 10x lens)