

# GBP4 antibody

Catalog Number: 110893

## Product name

GBP4 antibody

## Specificity

Human; other species not tested.

## Antibody description

GBP4 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, A431 cells. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 73 kDa

## Preparation

This antibody was obtained by immunization of GBP4 recombinant protein (Accession Number: NM\_052941). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110893(GBP4 antibody) at dilution of 1:600



IP Result of anti-GBP4 (IP:Catalog No:110893, 4ug; Detection:Catalog No:110893 1:500) with A431 cells lysate 2400ug.