

GDPD5 antibody



Catalog Number: 110931

Product name

GDPD5 antibody

Specificity

Human; other species not tested.

Antibody description

GDPD5 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, COLO 320 cells. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of GDPD5 recombinant protein (Accession Number: BC018771). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110931(GDPD5 Antibody) at dilution of 1:300