

GEFT-Specific antibody

Catalog Number: 110935

Product name

GEFT-Specific antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

GEFT-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HEK-293 cells, mouse heart tissue, mouse liver tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human pancreas tissue, human skeletal muscle tissue. Observed molecular weight by Western-blot: 64 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_182947). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations




mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110935(GEFT-Specific antibody) at dilution of



1:1500

IP Result of anti-GEFT-Specific (IP:Catalog No:110935, 4ug; Detection:Catalog No:110935 1:1500) with mouse brain tissue lysate 8000ug.




Immunohistochemistry of paraffin-embedded

GEFT-Specific antibody

Catalog Number: 110935

human pancreas slide using Catalog
No:110935(GEFT-Specific Antibody) at dilution of
1:100



Immunohistochemistry of paraffin-embedded
human pancreas slide using Catalog
No:110935(GEFT-Specific Antibody) at dilution of
1:100