

GFI1 antibody



Catalog Number: 110946

General Information

Product name:

GFI1 antibody

Immunogen:

Applications:

ELISA, WB

Specificity:

Human, Mouse, Rat; other species not tested.

Antibody description:

GFI1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HL-60 cells, PC-3 cells, Raji cells, RAW264.7 cells. Observed molecular weight by Western-blot: 62-66 kDa

Properties

Preparation:

This antibody was obtained by immunization of GFI1 recombinant protein (Accession Number: BC032751). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal


Ig Type: Rabbit IgG

Applications

Recommended Dilution:

WB: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110946(GFI1 antibody) at dilution of 1:500