

# GFPT2 antibody

Catalog Number: 110952

## Product name

GFPT2 antibody

## Specificity

Human; other species not tested.

## Antibody description

GFPT2 Rabbit Polyclonal antibody. Positive IP detected in HT-1080 cells. Positive WB detected in HT1080 cells. Observed molecular weight by Western-blot: 75-77 kDa

## Preparation

This antibody was obtained by immunization of GFPT2 recombinant protein (Accession Number: NM\_005110). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IP

## Dilutions


Recommended Dilution:

WB: N/AIP: 1:200-1:1000

## Validations



WB result of GFPT2 antibody (Catalog No:110952) with HT1080 cells.



IP Result of anti-GFPT2 (IP:Catalog No:110952, 4ug; Detection:Catalog No:110952 1:300) with HT-1080 cells lysate 2400ug.