

GFRA2 antibody

Catalog Number: 110953

Product name

GFRA2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GFRA2 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Positive IHC detected in human brain tissue, human placenta tissue. Observed molecular weight by Western-blot: 41-45kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001495). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110953(GFRA2 antibody) at dilution of 1:200



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110953(GFRA2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110953(GFRA2 Antibody) at dilution of 1:50 (under 40x lens)