

GGA1 antibody



Catalog Number: 110956

Product name

GGA1 antibody

Immunogen

[Human GGA1 Recombinant protein \(His tag\)](#)

Specificity

Human, Rat, Mouse; other species not tested.

Antibody description

GGA1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 70-72 kDa

Preparation

This antibody was obtained by immunization of GGA1 recombinant protein (Accession Number: XM_017028762). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110956(GGA1 Antibody) at dilution of 1:600