

GNAT1 antibody

Catalog Number: 110993

Product name

GNAT1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GNAT1 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000172). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IF: 1:50-1:500

Validations




WB results of GNAT1 antibody (Catalog No:110993) with WT mouse Eye and Prph2 (Rds) mutant mouse Eye (Negative control). Courtesy of Seongjin Seo, PhD, University of Iowa College of



Medicine.

mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:110993(GNAT1 antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:110993(GNAT1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)