

GNL1 antibody

Catalog Number: 111004

Product name

GNL1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

human; other species not tested.

Validations

Antibody description

GNL1 Rabbit Polyclonal antibody. Positive IHC detected in human heart. Positive IF detected in HeLa. Positive WB detected in MCF7 cells, HeLa cells, Jurkat cells, K-562 cells, Raji cells, U-937 cells. Observed molecular weight by Western-blot: 69-75kd




Preparation

This antibody was obtained by immunization of GNL1 recombinant protein (Accession Number: XM_005249015). Purification method: Antigen affinity purified.

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111004(GNL1 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human heart using Catalog No:111004(GNL1 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded human heart using Catalog No:111004(GNL1 antibody) at dilution of 1:100 (under 10x lens)

GNL1 antibody

Catalog Number: 111004

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using GNL1 antibody Catalog No:111004 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).