

GLT25D2 antibody



Catalog Number: 111038

General Information

Product name:

GLT25D2 antibody

Immunogen:

[Human GLT25D2 Recombinant protein \(GST tag\)](#)

Applications:

ELISA, WB

Specificity:

Human, Mouse; other species not tested.

Antibody description:

GLT25D2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, Neuro-2a cells. Observed molecular weight by Western-blot: 65 kDa

Properties

Preparation:

This antibody was obtained by immunization of GLT25D2 recombinant protein (Accession Number: BC035672). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:111038(GLT25D2 Antibody) at dilution of 1:600