

# GLYCTK antibody



Catalog Number: 111086

## Product name

GLYCTK antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GLYCTK Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer. Positive WB detected in human lung tissue, human liver tissue. Observed molecular weight by Western-blot: 55kd

## Preparation

This antibody was obtained by immunization of GLYCTK recombinant protein (Accession Number: BC021896). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111086(GLYCTK antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111086(GLYCTK antibody) at dilution of 1:50 (under 10x lens)