

# Glypican 5 antibody



Catalog Number: 111091

---

## Product name

Glypican 5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Glypican 5 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 64kd

## Preparation

This antibody was obtained by immunization of Glypican 5 recombinant protein (Accession Number: NM\_004466). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111091(GPC5 antibody) at dilution of 1:300