

GMPS antibody

Catalog Number: 111104

Product name

GMPS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GMPS Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HL-60 cells, Jurkat cells, mouse brain tissue, mouse testis tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of GMPS recombinant protein (Accession Number: NM_003875). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111104(GMPS antibody) at dilution of 1:800



IP Result of anti-GMPS (IP:Catalog No:111104, 4ug; Detection:Catalog No:111104 1:800) with HeLa cells lysate 1080ug.