

GPR120 antibody



Catalog Number: 111107

Product name

GPR120 antibody

Specificity

Human; other species not tested.

Antibody description

GPR120 Rabbit Polyclonal antibody. Positive IHC detected in human prostate hyperplasia tissue.

Preparation

This antibody was obtained by immunization of GPR120 recombinant protein (Accession Number: NM_001195755). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:111107(GPR120 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:111107(GPR120 Antibody) at dilution of 1:100 (under 40x lens).