

# GPR142 antibody

Catalog Number: 111114

---

## Product name

GPR142 antibody

## Specificity

Human; other species not tested.

## Antibody description

GPR142 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 51-55 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_181790). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:111114(GPR142 antibody) at dilution of 1:300