

# GPR171 antibody

Catalog Number: 111119

## Product name

GPR171 antibody

## Specificity

Human; other species not tested.

## Antibody description

GPR171 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue.

## Preparation

This antibody was obtained by immunization of GPR171 recombinant protein (Accession Number: NM\_013308). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111119(GPR171 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111119(GPR171 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).