

# GPR55 antibody



Catalog Number: 111131

## Product name

GPR55 antibody

## Specificity

Human; other species not tested.

## Antibody description

GPR55 Rabbit Polyclonal antibody. Positive IHC detected in human prostate hyperplasia tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005683).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:111131(GPR55 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:111131(GPR55 Antibody) at dilution of 1:200 (under 40x lens).