

# GPR65 antibody



Catalog Number: 111133

## Product name

GPR65 antibody

## Specificity

Human; other species not tested.

## Antibody description

GPR65 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue.

## Preparation

This antibody was obtained by immunization of GPR65 recombinant protein (Accession Number: NM\_003608). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111133(GPR65 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111133(GPR65 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).