

GPR83 antibody



Catalog Number: 111136

Product name

GPR83 antibody

Specificity

human; other species not tested.

Antibody description

GPR83 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 48-55kd

Preparation

This antibody was obtained by immunization of GPR83 recombinant protein (Accession Number: NM_016540). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111136(GPR83 antibody) at dilution of 1:1000