

GPRIN1 antibody

Catalog Number: 111141

Product name

GPRIN1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GPRIN1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 100-110 kDa

Preparation

This antibody was obtained by immunization of GPRIN1 recombinant protein (Accession Number: NM_052899). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111141(GRIN1 Antibody) at dilution of 1:600



IP Result of anti-GRIN1 (IP:Catalog No:111141, 3ug; Detection:Catalog No:111141 1:500) with mouse brain tissue lysate 3600ug.