

# GR repeat antibody



Catalog Number: 111153

## Product name

GR repeat antibody

## Specificity

Human; other species not tested.

## Antibody description

GR repeat Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein.

## Preparation

This antibody was obtained by immunization of Peptide. Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:111153(GR repeat Antibody) at dilution of 1:1000