

GR repeat antibody



Catalog Number: 111153

Product name

GR repeat antibody

Specificity

Human; other species not tested.

Antibody description

GR repeat Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of Peptide. Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:111153(GR repeat Antibody) at dilution of 1:1000