

GR repeat antibody



Catalog Number: 111153

General Information

Product name:

GR repeat antibody

Immunogen:

Applications:

ELISA, WB

Specificity:

Human; other species not tested.

Antibody description:

GR repeat Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein.

Properties

Preparation:

This antibody was obtained by immunization of Peptide. Purification method: Antigen Affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

Recommended Dilution:

WB: 1:500-1:5000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:111153(GR repeat Antibody) at dilution of 1:1000