

# GRAMD4 antibody



Catalog Number: 111156

## Product name

GRAMD4 antibody

## Specificity

Human; other species not tested.

## Antibody description

GRAMD4 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 66, 140 kDa

## Preparation

This antibody was obtained by immunization of GRAMD4 recombinant protein (Accession Number: XM\_047441209). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111156(GRAMD4 Antibody) at dilution of 1:300