

GSNOR,ADH5 antibody



Catalog Number: 111174

Product name

GSNOR,ADH5 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GSNOR,ADH5 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, human brain tissue, human heart tissue, human liver tissue, mouse brain tissue. Positive IHC detected in human prostate cancer. Observed molecular weight by Western-blot: 40kd

Preparation

This antibody was obtained by immunization of GSNOR,ADH5 recombinant protein (Accession Number: NM_000671). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:111174(GSNOR,ADH5 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:111174(GSNOR,ADH5 antibody) at dilution of 1:50 (under 40x lens)