

GST Tag antibody



Catalog Number: 111178

Product name

GST Tag antibody

Specificity

Tag

Antibody description

GST Tag Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: FN315687). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

Validations



Western blot of GST-tagged fusion protein with anti-GST-tag (Catalog No:117323) at various dilutions.