

# GSTA4 antibody

Catalog Number: 111182

## Product name

GSTA4 antibody

## Immunogen

[Human GSTA4 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

GSTA4 Rabbit Polyclonal antibody. Positive IP detected in SH-SY5Y cells. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 27 kDa

## Preparation

This antibody was obtained by immunization of GSTA4 recombinant protein (Accession Number: XM\_005249035). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111182(GSTA4 antibody) at dilution of 1:1200



IP Result of anti-GSTA4 (IP:Catalog No:111182, 4ug; Detection:Catalog No:111182 1:500) with SH-SY5Y cells lysate 1600ug.