

GRN1/G antibody



Catalog Number: 111216

Product name

GRN1/G antibody

Specificity

Human; other species not tested.

Antibody description

GRN1/G Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 64 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM_005257253). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:111216(GRN1/G Antibody) at dilution of 1:600