

# Growth hormone antibody



Catalog Number: 111218

## Product name

Growth hormone antibody

## Specificity

Human; other species not tested.

## Antibody description

Growth hormone Rabbit Polyclonal antibody.  
Positive WB detected in human placenta tissue.  
Observed molecular weight by Western-blot: 22kd

## Preparation

This antibody was obtained by immunization of Growth hormone recombinant protein (Accession Number: NM\_000515). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:111218(GH1 antibody) at dilution of 1:1000