

# GTF2H1 antibody



Catalog Number: 111232

---

## Product name

GTF2H1 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

GTF2H1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 62kd

## Preparation

This antibody was obtained by immunization of GTF2H1 recombinant protein (Accession Number: NM\_005316). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111232(GTF2H1 antibody) at dilution of 1:200