

# GTF3C4 antibody

Catalog Number: 111238

## Product name

GTF3C4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GTF3C4 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 92 kDa

## Preparation

This antibody was obtained by immunization of GTF3C4 recombinant protein (Accession Number: NM\_012204). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111238(GTF3C4



antibody) at dilution of 1:1600

IP Result of anti-GTF3C4 (IP:Catalog No:111238, 3ug; Detection:Catalog No:111238 1:800) with HeLa cells lysate 875ug.