

# GTPBP3 antibody



Catalog Number: 111243

## Product name

GTPBP3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GTPBP3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 45-50kd

## Preparation

This antibody was obtained by immunization of GTPBP3 recombinant protein (Accession Number: NM\_032620). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111243(GTPBP3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:111243(GTPBP3 antibody) at dilution of 1:100 (under 10x lens)