

GTPBP5 antibody

Catalog Number: 111245

Product name

GTPBP5 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Immunogen

[Human GTPBP5 Recombinant protein \(GST tag\)](#)


Validations

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GTPBP5 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 44 kDa




Preparation

This antibody was obtained by immunization of GTPBP5 recombinant protein (Accession Number: XM_011528776). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111245(GTPBP5 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded human testis using Catalog No:111245(GTPBP5 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:111245(GTPBP5

GTPBP5 antibody



Catalog Number: 111245

antibody) at dilution of 1:100 (under 40x lens)