

HCLS1 antibody



Catalog Number: 111283

Product name

HCLS1 antibody

Specificity

Human; other species not tested.

Antibody description

HCLS1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 75 kDa

Preparation

This antibody was obtained by immunization of HCLS1 recombinant protein (Accession Number: NM_005335). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111283(HCLS1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111283(HCLS1 Antibody) at dilution of 1:50 (under 40x lens)