

# Hemoglobin alpha antibody



Catalog Number: 111287

## Product name

Hemoglobin alpha antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Hemoglobin alpha Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in K-562 cells, mouse liver tissue, rat brain tissue, rat heart tissue. Observed molecular weight by Western-blot: 14-16 kDa, 27 kDa

## Preparation

This antibody was obtained by immunization of Hemoglobin alpha recombinant protein (Accession Number: NM\_000558). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:111287(HBA1 antibody) at dilution of 1:1000



IP Result of anti-HBA1 (IP:Catalog No:111287, 4ug; Detection:Catalog No:111287 1:500) with mouse liver tissue lysate 4000ug.