

Hemoglobin epsilon-Specific antibody



Catalog Number: 111289

Product name

Hemoglobin epsilon-Specific antibody

Specificity

Human; other species not tested.

Antibody description

Hemoglobin epsilon-Specific Rabbit Polyclonal antibody. Positive WB detected in recombinant protein. Observed molecular weight by Western-blot: 42kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005330). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



WB results of Catalog No:111289 (HBE1 antibody) at dilution of 1:2000.