

Hemopexin antibody



Catalog Number: 111290

Product name

Hemopexin antibody

Specificity

Human; other species not tested.

Antibody description

Hemopexin Rabbit Polyclonal antibody. Positive WB detected in human blood tissue, human placenta tissue. Observed molecular weight by Western-blot: 65-75kd

Preparation

This antibody was obtained by immunization of Hemopexin recombinant protein (Accession Number: BC005395). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:111290(HPX antibody) at dilution of 1:2000