

# HIP14; ZDHHC17 antibody



Catalog Number: 111353

## Product name

HIP14; ZDHHC17 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HIP14; ZDHHC17 Rabbit Polyclonal antibody.  
Positive IHC detected in human testis. Positive WB detected in Raji cells, human kidney tissue.  
Observed molecular weight by Western-blot: 70kd

## Preparation

This antibody was obtained by immunization of HIP14; ZDHHC17 recombinant protein (Accession Number: NM\_015336). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111353(HIP14; ZDHHC17 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human testis using Catalog No:111353(HIP14; ZDHHC17 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:111353(HIP14; ZDHHC17 antibody) at dilution of 1:50 (under 40x lens)