

# HIP1R antibody



Catalog Number: 111354

## Product name

HIP1R antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HIP1R Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HeLa cells, human liver tissue, human stomach tissue, mouse brain tissue, mouse heart tissue, mouse lung tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 120 kDa

## Preparation

This antibody was obtained by immunization of HIP1R recombinant protein (Accession Number: NM\_001303097). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:111354(HIP1R antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human placenta using Catalog No:111354(HIP1R antibody) at dilution of 1:100 (under 40x lens)

# HIP1R antibody

Catalog Number: 111354

---

IP Result of anti-HIP1R (IP:Catalog No:111354, 4ug; Detection:Catalog No:111354 1:1000) with mouse brain tissue lysate 4000ug.

