

HK1 antibody

Catalog Number: 111412

Product name

Polyclonal

HK1 antibody

Ig Type**Specificity**

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications**Antibody description**

ELISA, WB

HK1 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 103 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of HK1 recombinant protein (Accession Number: NM_000188). Purification method: Antigen affinity purified.

Validations**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111412(HK1 Antibody) at dilution of 1:300

Clonality