

HLA-A antibody

Catalog Number: 111417

Product name

HLA-A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

HLA-A Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human small intestine tissue. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 41-44kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002116). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111417(HLA-A antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:111417(HLA-A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:111417(HLA-A Antibody) at dilution of 1:50 (under 40x lens)